


Exp. 6**Isolation of Caffeine from Tea***Techniques*

Extraction, recrystallization

Equipment

Magnetic stirrer, vacuum filtration set up, rotavap

Chemicals

Sodium bicarbonate, sodium sulfate, dichloromethane, isopropanol, petrolether

Safety

Extraction with dichloromethane should be conducted under the hood

Procedure:

For 10 Tea-bags:

- Dissolve 20 g of sodium carbonate in 200 ml of water and fill it into the big flat beaker.

- Bring the solution to its boiling point, put ten of the tea bags into it and infuse them for ten minutes.
- Let it cool down for five minutes and continue cooling with ice-bath.
- Take the tea-bags off and press the remaining liquid out of them.
- Pour the liquid into a separatory funnel.
- Perform the extraction three times with 50 ml dichloromethane.
- Avoid too vigorous shaking, otherwise an emulsion is formed which is difficult to separate.
- Collect the extracts in an Erlenmeyer flask and dry it with some spatula-tips of sodium sulfate to get rid of water-traces.
- Filter the extract into a 250 ml round-bottom flask.
- Evaporate the solvent on the rotavap.
- Recrystallize the remaining solid from 5 ml isopropanol using the water bath.
- Vacuum filtrate the crystals formed.
- Wash it with some cold petrol ether.
- Let it dry on air, fill it into a vial and determine the weight

Tasks

1. Calculate the content in the tea used.
2. Determine the melting point